

POWERED BY STEMINA 🤍


Overview

GOALS:

1) Create metabolism-based diagnostic tests based on reproducible metabotypes 2) Examine the overlap of biochemical processes associated with <u>METHODS:</u> Quantitative LC-MS/MS measurement of plasma metabolites followed by statistical analysis for metabotypes associated with ASD using CAMP subjects **<u>RESULTS</u>: Reproducible metabotypes associated amine, energy, and purine metabolic</u>** processes were identified in 41% of CAMP ASD subjects


Introduction

- •The Children's Autism Metabolome Project (CAMP, ClinicalTrials.gov: NCT02548442) recruited 1100 children for metaboloimc analyses to identify altered metabolism associated with ASD
- •Utilized a metabotyping approach previously published (Biol. Psychiatry 2019; 85:345-354) that identified ASD subjects associated with branched-chain amino acid (BCAA) dysregulation
- •Metabotype-based stratification of ASD provides more biochemically homogenous subpopulations that offer the potential for tailored pharmacological, behavioral, and dietary interventions



Children's Autism Metabolome Project CAMP

- •8 clinical sites, research reliable ADOS-2, MSEL, and DSM-V
- •Large cohort of fasted subjects allows detection distinct ASD subpopulations using plasma
- Over 600 metadata items per subject from both subject and family medical histories
- Genetic samples available for analysis



Untargeted

Identification of Neurotransmitter-Associated and Other Metabotypes: Further Stratification of the Children's Autism

A. Smith¹, M. Ludwig¹, D. Brass¹, D. Sugden¹, L. Feuling¹, E. Donley¹, M. Natowicz², D. Amaral³, and R. Burrier¹ ¹Stemina Biomarker Discovery Inc., ²Cleveland Clinic, ³The M.I.N.D. Institute, University of California - Davis

Targeted Metabolomic Analysis Identified Potential Metabotypes

- Targeted non-quantitative analysis of 82 metabolites to identify biochemical areas associated with metabotypes
- Training set of CAMP ASD and TYP subjects evaluated

me	tab	oty	ype	? S


Metabotypes Identified in Targeted Metabolomic Analys							
Metabolite	SEN	SPEC	PPV	Biochemical Area			
α-Ketoglutarate	0.126	1.000	1.000	Energy Homeostasis			
Hypoxanthine	0.126	0.988	0.970	Purine Degradation			
Pyruvate	0.193	0.965	0.942	Energy Homeostasis			
Lactate	0.118	0.976	0.938	Energy Homeostasis			
GLY	0.098	0.976	0.926	Amine / Neuro			
Taurine	0.142	0.965	0.923	Purine Degradation			
ASP	0.094	0.976	0.923	Amine / Neuro			
ASN	0.098	0.988	0.962	Amine / Neuro			
Succinate	0.094	0.988	0.960	Energy Homeostasis			
PPV is reported for the CAMP population (precision)							

Metabotypes using Quantitative Assays

- Quantitative LC-MS/MS assays developed based on metabolomics results
- The CAMP samples were divided between a training set and an independent validation (test) set

Purine Degradation Related Metabotype

- •Metabolites measured: xanthine, hypoxanthine, inosine, uric acid, and taurine
- Hemolyzed samples with hemoglobin levels >200 mg/dl were excluded from analysis due to interference • A single reproducible metabotype was identified in 6.3% of CAMP ASD subjects



- Xanthine Oxidoreductase (XOR) is required for catabolism of purines
- •Elevated xanthine is correlated with taurine (p=.64), an amine which reduces XOR activity and is evaluated in MOCOS (molybdenum cofactor required by XOR) deficiency
- The data suggests that biomarkers of defective sulfite metabolism may provide a link to understanding the biology in a subset of children with ASD


References

Smith AM, King JJ, West PR, Ludwig MA, Donley ELR, Burrier RE, Amaral DG. Amino Acid Dysregulation Metabotypes: Potential Biomarkers for Diagnosis and Individualized Treatment for Subtypes of Autism Spectrum Disorder. Biol Psychiatry. 2019, 85(4)345:354; PMID: 30446206

Support and Funding


National Institute of Mental Health (NIMH) Grant 1R44MH107124-01, Nancy Lurie Marks Family Foundation, The Robert E. and Donna Landreth Family Fund


na Samples					
AII					
526					
213					
739					
71.2					
78.5					
59.6					
35.1 +/- 7.8					
32.5 +/- 8.8					
18 to 48					
62.9 +/- 17.7					
101.8 +/- 16.3					



- Metabolites measured: 32 amine containing and amino acid metabolites
- Test for amine related metabotypes not previously published (Smith et al. 2019)
- Five amine metabotypes identified 21.5% of CAMP ASD subjects
- •Altered neurotransmission, mitochondrial biology, nitrogen metabolism may be associated with metabotypes.







- Mitochondrial biology/energy metabolism
- The large ASD CAMP cohort allows detection of reproducible metabotypes with a prevalence > 5% • Biological processes associated metabotypes were identified in 41% of CAMP ASD subjects
- - Diverse amino acid pathways related to amino acid homeostasis
 - Purine catabolism
- Future Studies

 - Utilize additional computational approaches to identify metabolic changes predictive of ASD
 - Continue CAMP subject metadata for associations with metabotype positive populations
 - Launch clinical studies of dietary interventions within BCAA metabotype positive ASD subjects

α-Ketoglut

α-Ketoglι



Energy Homeostasis Related Metabotypes

•Metabolites measured: α -ketoglutarate, lactate, pyruvate, succinate, alanine and phenylalanine •Lactate, pyruvate, and alanine are commonly used to assess mitochondrial bioenergetic function

Ratio	SEN		SPEC		PPV		Diachamical Area	
Ratio	Train	Test	Train	Test	Train	Test	Biochemical Area	
etoglutarate / Lactate	0.074	0.134	0.981	0.981	0.905	0.946	Energy Homeostasis	
-Ketoglutarate / ALA	0.070	0.061	0.981	0.991	0.900	0.941	Energy Homeostasis	
Lactate / ALA	0.051	0.073	0.990	0.991	0.929	0.950	Energy Homeostasis	
Lactate / PHE	0.058	0.084	0.990	0.981	0.938	0.917	Energy Homeostasis	
Lactato - Pyruvate Alanine								

Glutamate α-Ketoglutarate

•22.3% of CAMP ASD identified by an energy related metabotype

Novel Amine Associated Metabotypes

Ratio	SEN		SPEC		PPV		Biochemical
Ralio	Train	Test	Train	Test	Train	Test	
SN / GLY	0.069	0.073	0.990	0.990	0.947	0.950	Amine / Neuro
LY / PHE	0.062	0.065	1.000	0.990	1.000	0.944	Amine / Neuro
IS / LEU	0.088	0.084	0.980	0.980	0.920	0.917	Amine / Neuro
/N / ORN	0.050	0.057	0.990	0.990	0.929	0.938	Amine / Neuro
/S / ORN	0.065	0.050	0.990	0.990	0.944	0.929	Amine / Neuro

Overlap of Metabotype Positive Subjects

- •41% of CAMP ASD subjects identified in the training and test sets • Metabotype positive subjects that share multiple metabotypes
- may have more complex metabolic phenotypes requiring further study
- Subjects identified by a single amine, purine, or energy process metabotype may have a specific metabolic dysregulation

Summary

• Develop additional metabolic panels based on ROS and gut microbiome associated metabolism